

How it works


The **LINK** skillfully maneuvers to offer optimal legroom and head room while it gracefully exits the vehicle and lowers the original seat to the desired height, for easy level transfers. If you prefer to stand when transferring, the **LINK** can assist you with its lift-up feature.

Features

- Easy to use
- Movements tailored to your needs
- Always use the original seat
- Personalized seating memory
- Unmatched legroom and comfort
- Forward and aft movement
- Adjustable tilt and lift of seat
- Seating behind the **LINK** unaffected
- Seatback reclines completely
- Retain heated seats in most models

Specifications

- Maximum capacity: 300 lbs
- Equipped with an emergency back-up system
- Power supplied by your vehicle
- Power is not drained while **LINK** is not in use
- Trim options: Black - Grey - Tan





LINK

Innovative vehicle access solution



Connecting you and your vehicle

The **LINK** is an innovation like no other. Its technology is unsurpassed. Keep the look, comfort and style of your vehicle. It is compact, reliable, safe and easy to use. Get your vehicle, your way.

Safety

We take our customers' safety seriously. We are dedicated to manufacturing quality mobility equipment. Our products undergo a strict certification process. Certification of our products involves different tests such as: crash testing, cycle testing and weight capacity testing. We are proud to offer a line of products that exceed the industry's standards.

Installation

The **LINK** is designed to fit in the driver or front passenger position as well as second row passenger position in a vast array of vehicles.

It is installed into the original seat mounts with no structural modifications to your vehicle. It can easily be removed to return the vehicle to its original condition without affecting the resale value.

Due to the diversity of mobility solutions available to you, please contact your local mobility dealer for a personalized evaluation and overview of options tailored to your needs.

Our products qualify as adaptive equipment and are eligible for reimbursement under the terms of Mobility programs offered by various vehicle manufacturers.

Dealer :

ADAPT SOLUTIONS
QUALITY MOBILITY EQUIPMENT

145 DAMASE-BRETON | SAINT-LAMBERT, QC, G0S 2W0

P 866 641 0419 | F 418 889 9838 | O 418 889 0419

www.adaptsolutions.com



creation www.mobilstudio.com

www.adaptsolutions.com